

LABCOLD
AUSTRALIA



Bariatric Mortuary Refrigerators

TGA
Health Safety
Regulation



Temperature

2°C to 4 °C



Berths

NMRB1 - 1

NMRB2 - 2



Dimensions



NMRB1 1140w x 2500d x 1370h

NMRB2 1140w x 2500d x 2010h

Features:

- Top Mounted System
- TGA Approved
- Frost free refrigeration
- Interior Lighting
- Cabinet Finish in Stainless Steel Interior/ Exterior
- Digital Control Panel with built in alarms
- Hi / Lo Read out



info@mediline.com.au

Specifications

Model	NMRB1	NMRB2
Doors	1	2
External Dimensions	1140w x 2500d x 1370h	1140w x 2500d x 2010h
Internal Dimensions	940w x 2350d x 850h	940w x 2350d x 1490h
Weight	175 Kg	285 Kg
Voltage	240 V	240 V
Watts	1872W	2050W
Refrigerant	R134A	R134A
Warranty	12 Months*	12 Months*
Trays	1	2
Trays size	890 x 2055	890 x 2055
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	Stainless Steel	Stainless Steel
Insulation	50 mm	50 mm
High/Low Temperature Alarm	✓	✓
Fan Forced Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door	✓	✓
Optional Equipment		
Port entry	✓	✓
Remote alarm	✓	✓

*Refer to warranty sheet

**Labcold reserves the right to alter the specifications without prior notice. Labcold is a Division of Mediline Australia Pty Ltd. **

7 Longford Ct Springvale VIC 3171
 P: +(61) 03 9562 3355
 info@mediline.com.au