



Blanket Warmers

Temperature

60°C



Capacity

BW5 80L
BW10 257L
BW21 365L



Dimensions

BW5 500w x 575d x 600h
BW10 850w x 575d x 900h
BW21 750w x 700d x 1800h



Features:

- Fan forced system
- TGA Approved
- Overheating alarm
- Double glazed doors
- Stainless int./ext. with chrome plated shelves
- Micro-processor control
- Temperature lock out
- Digital Display with Built In Alarms
- Hi / Lo Read out



Specifications

Model	BW5	BW10	BW21
Capacity	80 L	257 L	365 L
External Dimensions	500w x 575d x 600h	850w x 575d x 900h	750w x 700d x 1800h
Internal Dimensions	390w x 480d x 430h	740w x 485d x 590h	630w x 610d x 1590h
Weight	75. Kg	125 Kg	175 Kg
Voltage	240 V	240 V	240 V
Warranty	12 Months*	12 Months*	12 Months*
Shelves size	390 x 480	740 x 485	630 x 610
Number of Shelves	2	3	5
Digital Display	✓	✓	✓
Castors	X	✓	✓
Interior/Exterior Cabinet Finish	Stainless Steel	Stainless Steel	Stainless Steel
High/Low Temperature Alarm	✓	✓	✓
Glass Door	✓	✓	✓

*Refer to warranty sheet

**Labcold reserves the right to alter the specifications without prior notice. Labcold is a Division of Mediline Australia Pty Ltd. **

7 Longford Ct Springvale VIC 3171
P: +(61) 03 9562 3355
info@mediline.com.au