



Temperature



-18°C to -22 °C

Capacity



MLF115	115L
MLF400	350L
MLF600	570L

Dimensions




MLF115	600w x 600d x 830h
MLF400	600w x 620d x 1900h
MLF600	770w x 750d x 1900h

Labcold's MLF Range -20 °C of Spark Safe Medical & Laboratory Freezers provide the reliability that you need for your valuable items.

The static refrigerated shelves help keep even temperature throughout.

Features:

- Bottom Mounted System
- TGA Approved
- Static Refrigerated Shelves
- Audible Alarm
- Door Lock
- Digital Display with Built In Alarms
- Hi / Lo Read out



Specifications

	MLF115	MLF400	MLF600
Capacity	115 L	350 L	570 L
External Dimensions (including castors)	600w x 600d x 830h	600w x 620d x 1900h	770w x 750d x 1900h
Internal Dimensions	480w x 410d x 660h	480w x 450d x 1610h	580w x 660d x 1660h
Weight	55 Kg	75 Kg	95 Kg
Voltage	240V 50Hz	240V 50Hz	240V 50Hz
Watts	252w	720w	960w
Refrigerant	R134a	R134a	R134a
Warranty	12 Months *Refer to warranty sheet	12 Months *Refer to warranty sheet	12 Months *Refer to warranty sheet
Shelves size	500 x 500	475 x 410	575 x 620
Number of Shelves	2 Fixed Static	6 Fixed Static	6 Fixed Static
Digital Display	✓	✓	✓
Castors	X	Optional	Optional
Interior/Exterior Cabinet Finish	White Enamel	White Enamel	White Enamel
Insulation	75mm	60mm	60mm
High/Low Temperature Alarm	✓	✓	✓
Static Refrigeration System	✓	✓	✓
Door Lock	✓	✓	✓
Solid Door	✓	✓	✓
Independent data logger	✓	✓	✓
Optional Equipment			
Door ajar alarm	✓	-	-
Port entry	-	-	-
Voltage Free Contacts	✓	-	-
Temp alert	-	-	-
Chart recorder	-	-	-
Battery for backup controller	✓	-	-
Additional Shelves	-	-	-



**Labcold reserves the right to alter the specifications without prior notice. Labcold is a Division of Mediline Australia Pty Ltd. **