



## Temperature

2°C to 8°C

-18°C to -22°C



## Capacity

Refrigerator 386L

Freezer 212L



## Dimensions

NDT 870 x 750 x 2120



The NDT Spark Free contains a dual refrigeration system: a unit for refrigeration and another for freezer.

This locally manufactured cabinet offers quality 304 grade stainless steel interior and exterior as standard with the option of 316.

## Features:

- Top Mounted System
- TGA Approved
- Heated Cabinet facia to prevent condensation
- Auto Defrost in 2°C Section
- Manual Defrost in -20 °C Section
- Stainless Steel Interior / Exterior
- Digital Display with Built in Alarms
- Hi / Lo Read out



# Specifications

	Refrigerator Section	Freezer Section
Capacity	386L	212L
External Dimensions (including castors)	870w x 750d x 2120h	870w x 750d x 2120h
Internal Dimensions	770w x 725d x 885h	685w x 535d x 580h
Weight	125 Kg	125 Kg
Voltage	240 V	240 V
Power Requirement	10 amp GPO	10 amp GPO
Watts	960W	540W
Refrigerant	R134a	R404
Warranty	12 Months*	12 Months*
Temperature Range	2°C to 8°C	-18°C to -22°C
Shelves size	750 x 555	520 x 675
Number of Shelves	3	2
Digital Display	✓	✓
Castors	✓	-
Interior/Exterior Cabinet Finish	Stainless Steel	Stainless Steel
Insulation	50 mm	75 mm
High/Low Temperature Alarm	✓	✓
Static Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door	✓	✓
<b>Optional Equipment</b>		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temp alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓

\*Refer to warranty sheet

\*\*Labcold reserves the right to alter the specifications without prior notice. Labcold is a Division of Mediline Australia Pty Ltd. \*\*