

LABCOLD
AUSTRALIA



NHF - Spark Safe Freezer

Temperature

-18°C to -22°C



Capacity

NHF400 356L

NHF600 505L



Dimensions



NHF400 650w x 700d x 1855h

NHF600 800w x 750d x 1900h

The NHF Spark Safe Freezer for Medical and Laboratory is the smartest, economical and safety way to store your products.

This Freezer is proudly manufactured in Australia, providing the quality and reliability required.

Features:

- Top Mounted System
- TGA Approved
- Heated Cabinet facia to prevent condensation
- Auto Defrost
- Castors / Legs
- Stainless Steel Interior
- Digital Display with Built-in Alarms
- Adjustable Shelves
- Hi / Lo Read Out



info@mediline.com.au

Specifications

*Refer to warranty sheet

	NHF400	NHF600
Capacity	356 L	505 L
External Dimensions (without castors)	650w x 700d x 1855h	800w x 750d x 1900h
Internal Dimensions	470w x 535d	590w x 535d
Weight	95 Kg	115 Kg
Voltage	240 V	240 V
Power Requirement	10 amp GPO	10 amp GPO
Watts	780W	960W
Refrigerant	R404A	R404A
Warranty	12 Months*	12 Months*
Shelves size	450 x 515	570 x 515
Number of Shelves	5	5
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	Stainless Steel / White Colourbond	Stainless Steel / White Colourbond
Insulation	75 mm	75 mm
High/Low Temperature Alarm	✓	✓
Static Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door	✓	✓
Optional Equipment		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temp alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓