



Temperature



-18°C to -22 °C

Capacity



NHFTS400 350L

NHFTS600 570L

Dimensions



NHFTS400 650w x 750d x 1915h

NHFTS600 770w x 750d x 1965h

The NHFTS is a Top Mounted Freezer which has been designed to storage flammable materials safely. It meets with Zone 2 classification by being located 300mm above floor level.

Safe to use in research laboratories where Solvents, Flammable goods, Petroleum Bio Products, Chemicals, or other hazard items need to be stored.

Features :

- Top Mounted System
- TGA Approved
- Heated Cabinet facia to prevent condensation
- Auto Defrost
- Castors / Legs
- Stainless Steel Interior
- Digital Display with Built-in Alarms
- Hi / Lo Read Out



Specifications

	NHFTS400	NHFTS600
Capacity	350 L	570 L
External Dimensions (including castors)	650w x 750d x 1915h	770w x 750d x 1965h
Internal Dimensions	470w x 535d x 1405h	590w x 535d x 1455h
Weight	110 Kg	130 Kg
Voltage	240V 50Hz	240V 50Hz
Watts	720w	960w
Power Requirement	10amp GPO	10amp GPO
Power Supply	3 Pin plug	3 Pin plug
Refrigerant	R404a	R404a
Warranty	12 Months *Refer to warranty sheet	12 Months *Refer to warranty sheet
Shelves size	450 x 515	570 x 515
Number of Shelves	5 Adjustable	5 Adjustable
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	Int Stainless Steel / Ext. White Colourbond	Int Stainless Steel / Ext. White Colourbond
Insulation	75mm	75mm
High/Low Temperature Alarm	✓	✓
Static Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door	✓	✓
Optional Equipment		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temperature alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓
Additional shelves	✓	✓

*Refer to warranty sheet