

**LABCOLD**  
AUSTRALIA



# NHRS - Spark Safe Refrigerator

## Temperature

2°C to 8°C



## Capacity

NHRS400 350L

NHRS600 590L



## Dimensions

NHRS400 600w x 620d x 1850h

NHRS600 770w x 750d x 1900h



The NHRS Spark Safe Refrigerator for Medical and Laboratory is the smartest, economical and safety way to store your products.

This refrigerator is proudly manufactured in Australia, providing the quality and reliability required.

## Features:

- Bottom Mounted System
- TGA Approved
- E-Glass & Heated cabinet facia to prevent condensation up to 80% R.H.
- Glass or Solid Door
- Auto Water Disposal
- Castors
- Door Lock
- Digital Display with built in Alarms
- Hi / Lo Read out



[info@mediline.com.au](mailto:info@mediline.com.au)

# Specifications

## Optional equipment

Temperature Data Logger



Keypad

Model	NHRS400	NHRS600
Capacity	350 L	590 L
External Dimensions (without castors)	600w x 620d x 1850h	770w x 750d x 1900h
Internal Dimensions	500w x 470d x 1750h	670w x 600d x 1740h
Weight	75 Kg	95 Kg
Voltage	240 V	240 V
Power Requirement	10 amp GPO	10 amp GPO
Watts	480W	575W
Refrigerant	R134a	R134a
Warranty	12 Months*	12 Months*
Shelves size	480 x 450	650 x 580
Number of Shelves	5	5
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	White Colourbond	White Colourbond
Insulation	50 mm	50 mm
High/Low Temperature Alarm	✓	✓
Static Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door / Glass Door	✓	✓
<b>Optional Equipment</b>		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temp alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓

\*Refer to warranty sheet

\*\*Labcold reserves the right to alter the specifications without prior notice. Labcold is a Division of Mediline Australia Pty Ltd. \*\*

7 Longford Ct Springvale VIC 3171  
 P: +(61) 03 9562 3355  
 info@mediline.com.au