



NHRTS - Spark Safe Refrigerator

Temperature

2°C to 8°C



Capacity

NHRTS400 350L

NHRTS600 590L



Dimensions



NHRTS400 600w x 620d x 1915h

NHRTS600 770w x 750d x 1965h

The NHRTS Spark Safe Refrigerator for Medical and Laboratory is the smartest, economical and safety way to store your products.


This refrigerator is proudly manufactured in Australia, providing the quality and reliability required.

Features:

- Top Mounted System
- TGA Approved
- E-Glass & Heated cabinet facia to prevent condensation up to 80% R.H.
- Adjustable Shelves
- Door Lock
- Digital Display with Built In Alarms
- Hi / Lo Read out



Optional equipment



Model	NHRTS400	NHRTS600
Capacity	350 L	590 L
External Dimensions (including castors)	600w x 620d x 1915h	770w x 750d x 1965h
Internal Dimensions	500w x 470d x 1490h	670w x 600d x 1530h
Weight	95 Kg	115 Kg
Voltage	240V 50Hz	240V 50Hz
Watts	386w	480w
Power Requirement	10amp GPO	10amp GPO
Power Supply	3 Pin plug	3 Pin plug
Refrigerant	R134a	R134a
Warranty	12 Months*	12 Months*
Shelves size	480 x 450	650 x 580
Number of Shelves	5 Adjustable	5 Adjustable
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	White Colourbond	White Colourbond
Insulation	50mm	50mm
High/Low Temperature Alarm	✓	✓
Static Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door / Glass Door	✓	✓
Optional Equipment		
Door ajar alarm	-	-
Port entry	-	-
Remote alarm	-	-
Temperature alert	-	-
Independent data logger	-	-
Chart recorder	-	-
Battery backup for controller	-	-
Additional shelves	-	-

*Refer to warranty sheet

**Labcold reserves the right to alter the specifications without prior notice. Labcold is a Division of Mediline Australia Pty Ltd. **

7 Longford Ct Springvale VIC 3171
 P: +(61) 03 9562 3355
 info@mediline.com.au