



Temperature

2°C to 8°C



Capacity

NLMBS/2 820L

NLMBS/3 1171L



Dimensions

NLMBS/2 1200w x 700d x 1980h

NLMBS/3 1880w x 700d x 1980h



Select from our range of Australian made cabinets which provide the quality and reliability required to store your items.


The NLMBS Laboratory Spark Free is available in either glass or solid door as well as in different sizes to suit your needs. Choose from our standard range or ask for a make to order refrigerator to suit your needs.

Features:

- Bottom Mounted System
- TGA Approved
- E-Glass & Heated cabinet facia to prevent condensation up to 80% R.H.
- Glass or Solid Door
- Auto Water Disposal
- Castors
- Door Lock
- Digital Display with built in Alarms
- Hi / Lo Read out



Optional equipment



Model	NLMBS/2	NLMBS/3
Capacity	820 L	1171 L
External Dimensions (without castors)	1200w x 700d x 1980h	1880w x 700d x 1980h
Internal Dimensions	1100w x 600d x 1450h	1780w x 600d x 1450h
Weight	125 Kg	175 Kg
Voltage	240 V	240 V
Power Requirement	10 amp GPO	10 amp GPO
Watts	984W	1012W
Refrigerant	R134a	R134a
Warranty	12 Months*	12 Months*
Shelves size	1050 x 530	585 x 530
Number of Shelves	4	12
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	White Colourbond	White Colourbond
Insulation	50 mm	50 mm
High/Low Temperature Alarm	✓	✓
Static Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door / Glass Door	✓	✓
Optional Equipment		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temp alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓

*Refer to warranty sheet

**Labcold reserves the right to alter the specifications without prior notice. Labcold is a Division of Mediline Australia Pty Ltd. **

7 Longford Ct Springvale VIC 3171
 P: +(61) 03 9562 3355
 info@mediline.com.au