

LABCOLD
AUSTRALIA



NLMi Refrigerator Incubator

Temperature

10°C to 40°C



Capacity

NLMi 1000/2 1000L



Dimensions

NLMi 1000/2 1200w x 700d x 1980h




Features:

- Top mounted system
- TGA Approved
- Five (5) rows of PVC coated shelves behind each door
- Interior Lighting LED
- Castors
- Automatic water disposal
- Digital Display with Built in Alarms
- Hi / Lo Read out



info@mediline.com.au

Optional equipment



Model	NLMi 1000/2
Capacity	1000 L
External Dimensions (without castors)	1200w x 700d x 1980h
Internal Dimensions	1100w x 600d x 1555h
Weight	125 Kg
Voltage	240 V
Power Requirement	10 amp GPO
Watts	850W
Refrigerant	R134a
Warranty	12 Months *
Shelves size	1080 x 530
Number of Shelves	5
Digital Display	✓
Castors	✓
Interior/Exterior Cabinet Finish	White Colourbond
Insulation	50 mm
High/Low Temperature Alarm	✓
Fan Forced Refrigeration System	✓
Door Lock	✓
Solid Door / Glass Door	✓
Optional Equipment	
Door ajar alarm	✓
Port entry	✓
Remote alarm	✓
Temp alert	✓
Independent data logger	✓
Chart recorder	✓
Battery backup for controller	✓

*Refer to warranty sheet

**Labcold reserves the right to alter the specifications without prior notice. Labcold is a Division of Mediline Australia Pty Ltd. **

7 Longford Ct Springvale VIC 3171
 P: +(61) 03 9562 3355
 info@mediline.com.au