

LABCOLD
AUSTRALIA



NLMS - Laboratory Spark Safe

Temperature

2°C to 8°C



Capacity

NLMS700/1 700L

NLMS1000/2 1000L

NLMS1614/3 1614L



Dimensions

NLMS700/1 870w x 700d x 2100h

NLMS1000/2 1200w x 700d x 2100h

NLMS1614/3 1880w x 700d x 2100h



Select from our range of Australian made cabinets which provide the quality and reliability required to store your items.

The NLMS Laboratory Spark Free Refrigerator is available in either glass or solid door as well as in different sizes to suit your needs.

Features:

- Top Mounted System
- TGA Approved
- E-Glass & Heated cabinet facia to prevent condensation up to 80% R.H.
- Glass or Solid Door
- Auto Water Disposal
- Castors
- Door Lock
- Digital Display with built in Alarms
- Hi / Lo Read out



TGA
Health Safety
Regulation

info@mediline.com.au

Specifications

Optional equipment

Temperature Data Logger



Keypad

Model	NLMS 700/1	NLMS 1000/2	NLMS 1614 /3	NLMS 2170 /4
Capacity	700 L	1000 L	1614 L	2170 L
External Dimensions (without castors)	870w x 700d x 2100h	1200w x 700d x 2100h	1880w x 700d x 2100h	2400w x 700d x 1980h
Internal Dimensions	770w x 550d x 1555h	1100w x 550d x 1555h	1780w x 550d x 1520h	2300w x 550d x 1520h
Weight	95 Kg	125 Kg	175 Kg	195 Kg
Voltage	240V 50Hz	240V 50Hz	240V 50Hz	240V 50Hz
Power Requirement	10amp GPO	10amp GPO	10amp GPO	10amp GPO
Watts	850W	1000W	1128W	1970W
Refrigerant	R134a	R134a	R134a	R134a
Warranty	12 Months*	12 Months*	12 Months*	12 Months*
Shelves size	750 x 530	1080 x 530	565 x 530	575 x 530
Number of Shelves	5 Adjustable	5 Adjustable	15 Adjustable	20 Adjustable
Digital Display	✓	✓	✓	✓
Castors	✓	✓	✓	✓
Interior/Exterior Cabinet Finish	White Colourbond	White Colourbond	White Colourbond	White Colourbond
Insulation	50mm	50mm	50mm	50mm
High/Low Temperature Alarm	✓	✓	✓	✓
Static Refrigeration System	✓	✓	✓	✓
Door Lock	✓	✓	✓	✓
Solid Door / Glass Door	✓	✓	✓	✓
Optional Equipment				
Door ajar alarm	✓	✓	✓	✓
Port entry	✓	✓	✓	✓
Remote alarm	✓	✓	✓	✓
Temperature alert	✓	✓	✓	✓
Independent data logger	✓	✓	✓	✓
Chart recorder	✓	✓	✓	✓
Battery backup for controller	✓	✓	✓	✓
Additional shelves	✓	✓	✓	✓

*Refer to warranty sheet

**Labcold reserves the right to alter the specifications without prior notice. Labcold is a Division of Mediline Australia Pty Ltd. **

7 Longford Ct Springvale VIC 3171
 P: +(61) 03 9562 3355
 info@mediline.com.au