

LABCOLD
AUSTRALIA



NLMTS Top Mounted Spark Safe Refrigerator

Temperature

2°C to 8°C



Capacity

| | |
|--------------|-------|
| NLMTS 700/1 | 700L |
| NLMTS 1000/2 | 1000L |
| NLMTS 1614/3 | 1614L |
| NLMTS 2170/4 | 2170L |



Dimensions

| | |
|--------------|----------------------|
| NLMTS 700/1 | 870w x 700d x 1980h |
| NLMTS 1000/2 | 1200w x 700d x 1980h |
| NLMTS 1614/3 | 1880w x 700d x 1980h |
| NLMTS 2170/4 | 2400w x 700d x 1980h |

Features:

- Top Mounted System
- TGA Approved
- E-Glass & Heated cabinet facia to prevent condensation up to 80% R.H.
- Static Refrigeration System
- Glass or Solid Doors
- Castors
- Auto Defrost
- Digital display with built in alarms
- Hi / Lo Read out



info@mediline.com.au

Optional equipment



| Model | NLMTS 700/1 | NLMTS 1000/2 | NLMTS 1614/3 | NLMTS 2170/4 |
|---------------------------------------|---------------------|----------------------|----------------------|----------------------|
| Capacity | 700 L | 1000 L | 1614 L | 2170 L |
| External Dimensions (without castors) | 870w x 700d x 1980h | 1200w x 700d x 1980h | 1880w x 700d x 1980h | 2400w x 700d x 1980h |
| Internal Dimensions | 770w x 550d x 1555h | 1100w x 550d x 1555h | 1780w x 550d x 1520h | 2300w x 550d x 1520h |
| Weight | 95 Kg | 125 Kg | 175 Kg | 195 Kg |
| Voltage | 240 V | 240 V | 240 V | 240 V |
| Power Requirement | 10 amp GPO | 10 amp GPO | 10 amp GPO | 10 amp GPO |
| Watts | 850W | 1000W | 1128W | 1970 |
| Refrigerant | R134a | R134a | R134a | R134a |
| Warranty | 12 Months* | 12 Months* | 12 Months* | 12 Months* |
| Shelves size | 750 x530 | 1080 x 530 | 565 x 530 | 575 x 530 |
| Number of Shelves | 5 | 5 | 15 | 20 |
| Digital Display | ✓ | ✓ | ✓ | ✓ |
| Castors | ✓ | ✓ | ✓ | ✓ |
| Interior/Exterior Cabinet Finish | White Colourbond | White Colourbond | White Colourbond | White Colourbond |
| Insulation | 50 mm | 50 mm | 50 mm | 50 mm |
| High/Low Temperature Alarm | ✓ | ✓ | ✓ | ✓ |
| Static Refrigeration System | ✓ | ✓ | ✓ | ✓ |
| Door Lock | ✓ | ✓ | ✓ | ✓ |
| Solid Door / Glass Door | ✓ | ✓ | ✓ | ✓ |
| Optional Equipment | | | | |
| Door ajar alarm | ✓ | ✓ | ✓ | ✓ |
| Port entry | ✓ | ✓ | ✓ | ✓ |
| Remote alarm | ✓ | ✓ | ✓ | ✓ |
| Temp alert | ✓ | ✓ | ✓ | ✓ |
| Independent data logger | ✓ | ✓ | ✓ | ✓ |
| Chart recorder | ✓ | ✓ | ✓ | ✓ |
| Battery backup for controller | ✓ | ✓ | ✓ | ✓ |

*Refer to warranty sheet

**Labcold reserves the right to alter the specifications without prior notice. Labcold is a Division of Mediline Australia Pty Ltd. **

7 Longford Ct Springvale VIC 3171
 P: +(61) 03 9562 3355
 info@mediline.com.au