



Temperature

5°C to 80°C



Capacity

I9042 45L

I9082 80L



Dimensions


I9042 350w x 370d x 350h

I9082 400w x 405d x 500h



Features:

- Fan Forced System
- Auto controller of fan speed
- Shelf Check Function- easy to identify problems
- Over Temperature and temperature different alarms
- Programmable controller.
- Digital Display with Built In Alarms
- Hi / Lo Read out



Specifications

| Model | i9042 | i9082 |
|-------------------------------------------------|--------------------|--------------------|
| Capacity | 45 L | 80 L |
| Internal Dimensions | 350w x 370d x 350h | 400w x 405d x 500h |
| Voltage | 220V - 50Hz | 220V - 50Hz |
| Power Consumption | 250W | 300W |
| Timing Range | 1 ~ 9999min | 1 ~ 9999min |
| Number of Shelves | 2 | 2 |
| Digital Display | ✓ | ✓ |
| Temperature Range | RT +5~80°C | |
| Temperature Stability | ±0.3°C | |
| Ambient Temperature | +5~35°C | |
| Display Resolution | 0.1°C/±0.2°C | |
| Warranty | 12 Months* | |
| Optional Equipment | | |
| Intelligent Programmable Temperature Controller | ✓ | ✓ |
| Independent Limiting Alarm | ✓ | ✓ |
| R485 Connector | ✓ | ✓ |
| Printer | ✓ | ✓ |

*Refer to warranty sheet

**Labcold reserves the right to alter the specifications without prior notice. Labcold is a Division of Mediline Australia Pty Ltd. **

7 Longford Ct Springvale VIC 3171
 P: +(61) 03 9562 3355
 info@mediline.com.au